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1. Sepsis Fact Sheet. National Institute of General Medical Sciences, 2014. (Accessed Dec. 15, 2015, at https://www.nigms.nih.gov/Education/pages/factsheet_sepsis.aspx.)
2. GSA News and Info. Global Sepsis Alliance (GSA). (Accessed Dec. 15, 2015, at <http://globalsepsisalliance.com/gsa-news-and-info/>.)
3. Elixhauser A, Friedman B, Stranges E. Septicemia in U.S. Hospitals, 2009. HCUP Statistical Brief #122. Rockville, MD: Agency for Healthcare Research and Quality (AHRQ); 2011 October 2011.
4. Hall M, Williams S, DeFrances C, Golosinskiy A. Inpatient care for septicemia or sepsis: A challenge for patients and hospitals. NCHS data brief, no 62. Hyattsville, MD: National Center for Health Statistics; 2011.
5. B. A. Barnes, G.E. Behringer, F.C. Wheelock at all. Postoperative Sepsis: Trends and Factors Influencing Sepsis Over a 20-Year Period Reviewed in 20,000 Cases. /Annals of Surgery. Volume 154, Number 4, 1961. P. 585-598.
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC1465981/pdf/annsurg00858-0105.pdf>
6. T. Jakobsona, J. Karjagina, L. Vipp at all. Postoperative complications and mortality after major gastrointestinal surgery. /Medicina 50 (2014). P. 111-117.
<http://dx.doi.org/10.1016/j.medici.2014.06.002>
7. B.T. Bateman, U. Schmidt, M.F. Berman at all. Temporal Trends in the Epidemiology of Severe Postoperative Sepsis after Elective Surgery: A Large, Nationwide Sample. /Anesthesiology 4 2010, Vol.112, 917-925. doi:10.1097/ALN.0b013e3181cea3d0.
<http://anesthesiology.pubs.asahq.org/article.aspx?articleid=1932871>
8. Vogel TR, Dombrovskiy VY, Lowry SF. Trends in postoperative sepsis: are we improving outcomes? / Surg Infect (Larchmt). 2009 Feb;10(1):71-8. doi: 10.1089/sur.2008.046.
<https://www.ncbi.nlm.nih.gov/pubmed/19298170>
9. T. R. Vogel, V. Y. Dombrovskiy, J. L. Carson at all. Postoperative sepsis in the United States./ Ann Surg. 2010 Dec; 252(6): 1065–1071. doi: 10.1097/SLA.0b013e3181dcf36e.
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2951484/>
10. Sepsis contributes to as many as half of all hospital deaths in U.S., study says. / Health System. University of Michigan. <http://www.uofmhealth.org/news/archive/201405/sepsis-contributes-many-half-all-hospital-deaths-us-study>
11. Celeste M. Torio, Brian J. Moore. National Inpatient Hospital Costs: The Most Expensive Conditions by Payer, 2013. STATISTICAL BRIEF #204, May 2016.
https://www.hcup-us.ahrq.gov/reports/statbriefs/sb204-Most-Expensive-Hospital-Conditions.jsp?utm_source=AHRQ&utm_medium=AHRQSTAT&utm_content=Content&utm_term=HCU P&utm_campaign=AHRQ_SB_204_2016.
12. Vincent Liu, MD, MS Gabriel J. Escobar at all. Hospital Deaths in Patients with Sepsis From 2 Independent Cohorts. JAMA July 2, 2014 Volume 312, Number 1, p.90-92. <http://jamanetwork.com/journals/jama/fullarticle/1873131#.U3kivi3tbAw.twitter>.
13. CDC. Making Health Care Safer - Think sepsis. Time matters. CDC Vital Signs 2016.
14. Karidis A. Sepsis is a medical emergency, CDC says. It can be stopped if caught in time. The Washington Post 2016 August 23, 2016.
15. Davis CP. What are the signs or symptoms of sepsis (blood poisoning)? In: MedicineNet, Inc.; 2015.
16. Jin M, Khan AI. Procalcitonin: Uses in the Clinical Laboratory for the Diagnosis of Sepsis. Laboratory Medicine 2010;41:173-7.
17. Gelfand B, Filimonov M, Hyles T. Procalcitonin: a new laboratory marker of sepsis and septic complications in surgery. Journal of Intensive Care 2003;1.

18. Pierrakos C, Vincent JL. Sepsis biomarkers: A review. Crit Care. 2010;14(1):R15. PMC2875530
19. T. Davydova, I. Matveeva, A. Gritsai, N. Pogosyan, L. Mamedova, V. Kuznetsov. Clinical and diagnostic value of the functional and conformational properties of serum albumin in patients with ovarian cancer (a pilot study). Women Reproductive System Tumors, 2014, 2, p.59-61 (*in Russian*).
20. Kalachik A., Ugolev I., Zabello T., Sadosky D., Oganova E., Murauski U. Electron Spin Resonance Spectroscopy Evaluation of the Serum Albumin Function Quality in Kidney Transplant Recipients. Proceedings of the National Academy of Sciences of Belarus, Series of Medical Sciences, 2014, Volume 4, p. 72-77 (*in Russian*).
21. A. Kalachyk, I. Ugolev, T. Zabello, E. Oganova, V. Muravsky. Electron Spin Resonance Spectroscopy of Serum Albumin is a Novel Test for Diagnosis of the Kidney Transplant Acute Rejection. Clin Chem Lab Med 2014; 52, Special Suppl, pp S1 – S1760, June 2014, P. 1068.
22. P. Vorobiev, L. Bezmelnitsyna, L. Krasnova, M. Holownia, I. Matveeva, T. Davydova, I. Nekhaev, A. Sytov, M. Nersesyan. Clinical and economical analysis of efficiency of ATA-test for early diagnosis of postoperative purulent-septic complications. Health care standardization problems, 2014, Volume 3-4, pp. 28-36 (*in Russian*).
23. Vorobiev P., Bezmelnitsyna L., Krasnova L., Muravsky V., Davidova T., Sytov A. Economical Justification of the Innovation Method of Laboratory Diagnostic of the Structural and Functional Changes of Serum Albumin in Septic Complications (Ata-Test). Value in Health, 2013, Volume 16, Issue 7, Page A364.
24. V. Sidorenko, V. Muravsky, E. Shimanovich, L. Zenko. Deterioration of serum albumin transport function in pregnant women with gestosis. Medical Journal (BSMU, Belarus), 2013, No 1, 149-150 (*in Russian*)
25. Liu et al. Changes in the electron paramagnetic resonance spectra of albumin-associated spin-labeled stearic acid as a diagnostic parameter of colorectal cancer. World Journal of Surgical Oncology, 2013, 11:223.
26. Davidova T., Matveeva I., Citov A., Nekhaev I., Sviridova S., Zubrikhina G., Jujginova O., Muravsky V., Gurachevskaya T. Diagnostics and monitoring of a development of purulo-septic complications in early postsurgical period in cancer patients using EPR method. International Journal of Applied and Fundamental Research. 2011, No 5, 94-95 (*in Russian*).
27. R. Jalan, K. Schnurr, R.P Mookerjee, S. Sen, L. Cheshire, S. Hodges, V. Muravsky, R. Williams, G. Matthes, N. A Davies. Alterations in the functional capacity of albumin in patients with decompensated cirrhosis is associated with increased mortality. Hepatology. 2009 Aug; 50(2); 555-64.
28. V. Muravsky, T. Gurachevskaya, S. Berezenko, K. Schnurr, A. Gurachevsky. Fatty acid binding sites of human and bovine albumins: Differences observed by spin probe ESR. Spectrochimica Acta Part A: Molecular and Biomolecular Spectroscopy, 74 (2009).
29. A. Gurachevsky, S. C. Kazmierczak, A. Jörres and V. Muravsky. Application of Spin Label EPR in the diagnosis and prognosis of cancer and sepsis. Clin Chem Lab Med 2008; 46(9); 1203-1210
30. V. Muravsky, A. Gurachevsky, G. Matthes. Disease-specific albumin patterns defined by electron spin resonance. Tumor Biol. 28 (suppl. 1), 19 (2007). ISSN 1010-4283; 1423-0380
31. Gurachevsky A, Muravskaya E, Gurachevskaya T, Smirnova L, Muravsky V. Cancer-Associated Alteration in Fatty Acid Binding to Albumin Studied by Spin-Label Electron Spin Resonance. Cancer Invest. 2007. 25:378-383. ISSN: 0735-7907
32. A. Gurachevsky, E. Shimanovitch, T. Gurachevskaya, V. Muravsky. Intra-albumin migration of bound fatty acid probed by spin label ESR. Biochem Biophys Res Commun. 2007 Sep 7; 360(4):852-856
33. S. C. Kazmierczak, A. Gurachevsky, G. Matthes, and V. Muravsky. Electron Spin Resonance Spectroscopy of Serum Albumin: A Novel New Test for Cancer Diagnosis and Monitoring. Clinical Chemistry 52:11, 2129–2134 (2006)

34. P.Seidel, A.Gurachevsky, V.Muravsky, K.Schnurr, G.Seibt, G.Matthes. Recognition of malignant processes with neural nets from ESR spectra of serum albumin. Z. Med. Phys. 15 (2005) 265–272
35. K. Schnurr, V. Muravsky, G. Seibt, G. Matthes. The analysis of functional characteristics of serum albumin – a novel possibility of early cancer diagnosis. In the book “Methods, possibilities, and perspectives of pre-symptomatic tumor diagnostics”, Pabst Science Publishers, 2005, p. 225-234. ISBN 1-59326-077-6 (USA)
36. Muravskaya E., Lapko A., Muravskii V. Modification of transport function of plasma albumin during atherosclerosis and diabetes mellitus. Bull Exp Biol Med. 2003 May; 135(5):433-5 (*in Russian*)
37. G.Matthes, G.Seibt, V.Muravsky, G.Hersmann, G.Dornheim. Albumin transport analysis of different collected and processed plasma products, by electron spin resonance spectroscopy. Transfus Apheresis Science, 2002, 27, p. 129-135
38. E.Melnichenko, A.Kuschner, A.Milutin, A.Mashevsky, V.Muravsky. Evaluation by the enamel of temporary teeth of the ionizing radiation doses adsorbed by the children population of Republic of Belarus. Stomatologiya (Moscow) 2002, 4, p. 54-56 (*in Russian*)
39. A.Mashevsky, V.Prokchorova, V.Muravsky, S.Linev, I.Ugolev. Technology of the tooth enamel ESR-dosimetry for the individual adsorbed dose estimation. Ukrainian Radiological Journal, 1999, № 1, p. 49-50 (*in Ukrainian*)
40. V.Prokchorova, A.Mashevsky, T.Zirus, S.Lappo, L.Paschko, R.Smoliakova, L.Derjavez, N.Jurevich, V.Muravsky. The use of laboratory testing methods for breast cancer diagnostics. In the book "Actual problems of oncology and clinical radiology". Research Institute of oncology and clinical radiology, Minsk, 1998, p. 77-81 (*in Russian*)
41. R.Smoliakova, A.Mashevsky, V.Prokchorova, T.Zirus, L.Paschko, A.Milutin, V.Muravsky. Investigation of structural and functional status of albumin in cancer patients by means of ESR-spectroscopy method. Advances of medical science in Belarus, 1997, 2, p. 19-20 (*in Russian*)
42. A.Mashevsky, V.Muravsky, E.Korotkevich, S.Linev, Y.Darvin, I.Ugolev, E.Javrid, V.Prokchorova, J.Zmuschko. Methodical and metrological aspects of gamma-irradiation adsorbed dose-measuring by tooth enamel. Medical physics, 1996, № 3, p. 46-50 (*in Russian*)
43. V.Prokchorova, A.Mashevsky, V.Muravsky, A.Milutin, T.Zirus, S.Linev, R.Smoliakova, L.Paschko. Diagnostic potential of ESR-spectrometric determination of conformational-dynamical characteristics of albumin for breast cancer. In the book “Actual problems of oncology and clinical radiology”. Research Institute of oncology and clinical radiology, Minsk, 1996, p. 47-55 (*in Russian*)
44. V.Muravsky A.Mashevsky, N.Krisenko, E.Javrid, S.Linev, I.Ugolev, Y.Darvin, V.Prokchorova. Parameters of the ESR spectrum spin model of tooth enamel and the use of it in ESR-dosimetry. Ukrainian Radiological Journal. 1995, № 4, p. 353-357 (*in Ukrainian*)
45. S.Linev, Y.Darvin, A.Mashevsky, V.Muravsky. Estimation of gamma-radiation absorbed doses by tooth enamel using spin model of ESR spectrum. Medical physics, 1995, №2, p. 14. (*in Russian*)
46. E.Javrid, A.Mashevsky, V.Muravsky, V.Prokchorova, E.Korotkevich, A.Milutin. A comparative study of conformation properties of human serum albumin in lung cancer by spin probe technique. Annals of Oncology, 1994, vol.5 (suppl.8), p. 9
47. A.Mashevsky, V.Muravsky, V.Prokchorova, A.Milutin. Comparative investigation of long chain fatty acid binding sites of normal serum albumin with the one from lung cancer patients carried out by the spin probe method. In the book “Blood serum albumin in clinical medicine”. Serie “Fluorescence methods of investigation and clinical diagnostics”, Edited by J.Grizonov and E.Dobretsov, Moscow, “IRUS”, 1994, p. 176-181 (*in Russian*)
48. A.Mashevsky, V.Muravsky, A.Milutin, E.Javrid, E.Korotkevich, V.Prokchorova, T.Zirus, S.Linev. Dynamics of ESR spectra of the spin-labeled blood serum under neoplastic process. Ukrainian Radiological Journal. 1994, № 4, p. 268-270 (*in Ukrainian*)

49. A.Mashevsky, V.Muravsky. The algorithms of bio-molecules parameters estimation for the development of clinical diagnostics procedures on the base of spin probe ESR method. In the book "Actual problems of oncology and clinical radiology", Research Institute of Oncology and Clinical Radiology, Belarus Ministry of Health, Minsk, 1993. P. 31-33 (*in Russian*)
50. O. Liang, K. Schnurr, K. Waterstradt. Laboratory diagnosis of active malignant proliferation by assessment of changes to transport parameters of serum albumin using EPR spectroscopy. Laboratory service. 2015;4(3): 44-56 (*in Russian*)
51. D.Nasonova, T.Ivanetz. "Deterioration of serum albumin transport function in preeclampsia"// Reference Book of the Head of Clinical Diagnostic Laboratory. June 2015, No 6, pp. 35-40 (*in Russian*)

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EN 55011-2012	Industrial, Scientific And Medical Equipment - Radio-Frequency Disturbance Characteristics - Limits And Methods Of Measurement
IEC 60065-2004	Audio, Video And Similar Electronic Apparatus - Safety Requirements
IEC 60601-1-2-2006	Medical electrical equipment - General requirements for basic safety and essential performance
IEC 61000-3-2-2006	Electromagnetic compatibility (EMC) - Limits - Limits for harmonic current emissions (equipment input current ≤ 16 A per phase)
IEC 61000-3-3-2005	Electromagnetic compatibility (EMC) - Limits - Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection.
IEC 61000-4-2-2006	Electromagnetic compatibility (EMC) - Testing and measurement techniques - Electrostatic discharge immunity test
IEC 61000-4-3-2009	Electromagnetic compatibility (EMC) - Testing and measurement techniques - Radiated, radio-frequency, electromagnetic field immunity test
IEC 61000-4-4-2006	Electromagnetic compatibility (EMC) - Testing and measurement techniques - Electrical fast transient/burst immunity test
IEC 61000-4-5-2006	Electromagnetic compatibility (EMC) - Testing and measurement techniques - Surge immunity test
IEC 61000-4-6-2011	Electromagnetic compatibility (EMC) - Testing and measurement techniques - Immunity to conducted disturbances, induced by radio-frequency fields
IEC 61000-4-8-2006	Electromagnetic compatibility (EMC) - Testing and measurement techniques - Power frequency magnetic field immunity test
IEC 61000-4-11-2006	Electromagnetic compatibility (EMC) - Testing and measurement techniques - Voltage dips, short interruptions and voltage variations immunity tests
EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control, and laboratory use - General requirements
EN 61010-2-081:2002	Safety requirements for electrical equipment for measurement, control, and laboratory use - Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes
IEC 61326-1:1997	Electrical equipment for measurement, control and laboratory use - EMC requirements - General requirements